

Product Information

Sperm Lysis Buffer, 50 mL

C/No. K.BUF-9106-50ml
Packaging 50 mL
Storage Room Temperature
(21 °C – 25 °C)

Apical Scientific Sdn Bhd
No. 17, Jalan BS 7/1C,
Taman Perindustrian Bukit
Serdang, 43300
Seri Kembangan, Selangor



Page 1 of 3

Description

Sperm Lysis Buffer is a solution formulated to lyse sperm cells and degrade proteins at the presence of Proteinase K and DTT (not provided – to be added fresh). This buffer is primarily used for genomic DNA extraction for downstream analysis.

Applications

- Genomic DNA extraction

Page 2 of 3

Quality Control

Assay:

DNA purity

Specification:

Absorbance with the ratio 260/280 between 1.7 – 2.0 was obtained. DNA band was observed.

Page 3 of 3