

Product Information

Stool DNA Collection Kit, 10 mL

C/No.	K.BUF-9107
Packaging	10 mL (1 prep) Available in (2 preps/ 10 preps/ 30 preps)
Storage	Room Temperature (21 °C – 25 °C)

Apical Scientific Sdn Bhd
No. 17, Jalan BS 7/1C,
Taman Perindustrian Bukit
Serdang, 43300
Seri Kembangan, Selangor



Page 1 of 3

Description

The Stool DNA Collection Kit is intended for the collection, preservation, and storage of stool samples for subsequent DNA extraction. It contains a non-lysis preservation solution that stabilizes nucleic acids and prevents degradation during sample transport or storage. This ensures reliable DNA recovery for downstream molecular applications.

Applications

- Genomic DNA extraction.

Page 2 of 3

Quality Control

<u>Assay:</u>	<u>Specification:</u>
Appearance	Clear, colorless Liquid
pH (neat; 25 °C)	11.2 – 11.50
Conductivity	7.7 – 8.0 mS/cm

Page 3 of 3