

Product Information

Lithium Chloride Solution, 8 M

C/No. K.RGT-9102-10ml
Packaging 10 mL
Storage Room temperature
(21 °C – 25 °C)

Apical Scientific Sdn Bhd
No. 17, Jalan BS 7/1C,
Taman Perindustrian Bukit
Serdang, 43300
Seri Kembangan, Selangor



Page 1 of 3

Description

Lithium Chloride Solution is a ready-to-use and nuclease-free solution for precipitating RNA. It is a fast and efficient method for purifying RNA from impurities such as DNA, protein or carbohydrates.

Applications

- RNA precipitation

Concentration

8 mol/L

Page 2 of 3

Quality Control

Assay:

DNase / RNase Activity

Specification:

None detected

Page 3 of 3