


Product Information

3M Sodium Acetate, pH5.2



C/No.	K.RGT-9106-1.5ml
Concentration	3M
Packaging	1.5 mL
Storage	Room Temperature (21 °C – 25 °C)

Apical Scientific Sdn Bhd
No. 17, Jalan BS 7/1C,
Taman Perindustrian Bukit
Serdang, 43300
Seri Kembangan, Selangor



Page 1 of 3


Description

3M Sodium Acetate is a 0.2 µm membrane filtered, nuclease-free solution. The solution is commonly used in nucleic acid extraction and purification workflows.

Grade

Biotechnology

Applications

- DNA precipitation
 - RNA precipitation
- 

Page 2 of 3

Quality Control

Assay:

DNase / RNase Activity

Specification:

None detected



Page 3 of 3