

Product Specification Sheet



Product Code: CUL-1101
Identifier: DMEM, high glucose, with L-Alanyl-L-Glutamine and Sodium Pyruvate
Grade: Cell Culture Media
Appearance: Clear Liquid
Molecular Formula: Not Applicable
Molecular Weight: Not Applicable
Description: Not Applicable
Storage Condition: Store between 2°C - 8°C. Protect from light.
Additional Info: Sterilized by 0.2µm membrane filtration.
Shelf Life: 12 Months

Assay:	Specification:
Appearance	Clear Liquid
pH (neat, 25°C)	7.00 - 7.50
Endotoxin activity (EU/ml)	< 0.25
Microbial Activity Test (5 days)	PASS
Osmolality (mOsm/kg)	320 - 360

Axil Scientific Pte Ltd

2 Tukang Innovation Grove, #06-01 JTC MedTech Hub, Singapore 618305
Tel: +65-6775 7318 Email: custcare@axilscientific.com
Website: www.axilscientific.com