

Product Specification Sheet



Product Code : BUF-1170

Identifier : Water, DEPC treated

Grade : Biotechnology

Appearance : Clear, colorless liquid

Molecular Formula : H₂O

Molecular Weight : 18.02 g/mol

Product Description : Water incubated with Diethylpyrocarbonate (DEPC) and autoclaved to remove DEPC.

Storage Condition : Store at room temperature.

<u>Assay:</u>	<u>Specification:</u>
Appearance	Clear Liquid
DNase/ RNase activity	None Detected
Protease Test	None Detected

Axil Scientific Pte Ltd

2 Tukang Innovation Grove #06-01, JTC MedTech Hub, Singapore 618305

Tel: +65-6775 7318 Email: info@axilscientific.com

Website: www.axilscientific.com