

Product Specification Sheet



Product Code : BUF-4100

Identifier : Sterile Deionized Water

Grade : Biotechnology

Appearance : Clear, colorless liquid

Molecular Formula : H₂O

Molecular Weight : 18.02g/mol

Product Description : Sterilized by 0.1µm membrane filtration. Suitable for molecular biology applications.

Storage Condition : Store at room temperature.

<u>Assay:</u>	<u>Specification:</u>
Appearance	Clear Liquid
DNase/ RNase activity	None Detected
Protease Test	None Detected

Axil Scientific Pte Ltd

2 Tukang Innovation Grove #06-01, JTC MedTech Hub, Singapore 618305
Tel: +65-6775 7318 Email: info@axilscientific.com
Website: www.axilscientific.com